

Elite Mini PAP Pen

- Catalog No.:** K 042
- Intended Use:** To provide a hydrophobic barrier when a circle is drawn around a specimen on the glass slide.
- Introduction:** Elite PAP pen provides a thin, film like water repellent barrier when a circle is drawn around the specimen on the glass slide. This barrier creates the proper surface tension to hold antibody solution within the target area on the slide. This avoids the repeated wipeout around the specimen during the staining procedure. Since the reagent is guarded by the hydrophobic circle around the tissue, it minimizes the undesirable dilution of the reagents with the wash buffer. This pen can be used for PAP, ABC, immunofluorescence method, frozen sections, and in situ hybridization procedures. The Mini PAP pen is especially useful when working with multiple sections on a single slide because its tip is smaller than PAP pen.
- Properties:**
1. Heat resistance up to 120°C.
 2. Ready to be used directly on wet slides.
 3. Applicable to microwave staining and in situ hybridization.
 4. Circle made by Elite PAP pen is insoluble in alcohol and acetone.
 5. Withstands rehydration and dehydration steps in alcohol.
 6. If desired, circle can be removed by xylene after the staining procedure is complete.
- Preparation before use:** Before use, hold the pen point up and press it you're your fingertip to take blocked air out from the pen barrel to prevent overflow.
- Procedure:** To start, press the tip of the pen straight down until fluid starts soaking the tip and then draw a circle around the specimen.
- Storage:** Always keep the pen tightly capped and store horizontally. Press down on the tip when pen starts to dry out.

IVD: For In Vitro Diagnostic Use

DBS will not be held responsible for patent infringement or other violation that may occur with the use of our product

DBS

1020 Serpentine Lane, # 114, Pleasanton, CA 94566 Tel: 925 484 3350, Fax: 925 484 3390

Website: www.dbiosys.com e-mail: customersupport@dbiosys.com