

Monoclonal Mouse Antibody to ZAP-70

Catalog No.:	PDM 149
Intended Use:	This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy. Clinical interpretation of staining results should be accompanied by histological studies with proper controls. Patients' clinical histories and other relevant diagnostic tests should be utilized by a qualified person(s) when evaluating and interpreting results.
Immuogen:	GST-fusion to tandem SH2 domains of human ZAP-70 corresponding to residues 1-254.
Clone:	2F3.2
Isotype:	IgG2a
Format:	This antibody has been pretitered and quality controlled to work on formalin-fixed paraffin-embedded and acetone fixed cryostat tissue sections. No further titration is required.
Staining Protocol:	We suggest an incubation period of 60 minutes at room temperature. Optimal incubation conditions should be determined by the user based upon the fixation conditions and staining system employed. <u>Formalin fixed paraffin embedded tissue sections require high temperature antigen unmasking with 10mM citrate buffer, pH 6.0 prior to immunostaining.</u>
Specificity:	Mutations in ZAP-70 gene results in a form of Severe Combined Immunodeficiency Syndrome (SCID) in humans. ZAP-70 expression also defines a subset of Chronic Lymphocytic Leukemia (CLL) in patients with unmutated Ig gene and poor clinical course. Recent studies suggest that protein levels of ZAP-70 are elevated in B-cells of CLL patients with non-mutant heavy-chain variable region (IgVH) but not those with the mutant regions. Recent evidence suggests that ZAP-70 could be an excellent prognostic biomarker with high levels of the proteins indicating a poor prognosis.
Positive Control:	Tonsil
Cellular Localization:	Cytoplasmic
Storage:	Store at 2-8°C. Do not use beyond the expiration date stated on the label.
References:	Crespo et al. Engl J Med 348 : 1764, 2003 Silvin et al. J Biol Chem 276:21450, 2001 Wiester et al. Blood 101: 4944, 2003.

IVD: For In Vitro Diagnostic Use

DBS will not be held responsible for patent infringement or other violation that may occur with the use of our product

DBS

1020 Serpentine Lane, # 114, Pleasanton, CA 94566 Tel: 925 484 3350, Fax: 925 484 3390

Website: www.dbiosys.com e-mail: customersupport@dbiosys.com