

Rabbit Polyclonal Antibody to *Helicobacter pylori*

Catalog No.:	PDR 053
Applications:	Analyte Specific Reagent. Analytical and Performance characteristics are not established. Optimal working conditions must be determined by the end user.
Immunogen:	Total lysate of <i>H. pylori</i> was injected in the rabbit.
Host:	Rabbit
Specificity:	This antibody reacts with whole <i>H. pylori</i> . It stains bacteria present on the surface of the epithelium or in the cytoplasm of the infected epithelial cells.
Positive Control:	Infected stomach
Cellular Localization:	Cytoplasmic
Storage:	Store at 2-8°C. Do not use beyond the expiration date stated on the label.
References:	i) Anderson et al. Apmis 96: 559, 1988.

ASR: Analyte Specific Reagent

DBS will not be held responsible for patent infringement or other violation that may occur with the use of our product

DBS

1020 Serpentine Lane, # 114, Pleasanton, CA 94566 Tel: 925 484 3350, Fax: 925 484 3390

Website: www.dbiosys.com e-mail: customersupport@dbiosys.com