

## Rabbit Monoclonal Antibody to Human c-erbB-2 Protein

Catalog No.:	RMAB 008
Applications:	Analyte Specific Reagent. Analytical and Performance characteristics are not established. Optimal working conditions must be determined by the end user.
Clone:	SP3
Immunogen:	Recombinant protein encoding extracellular domain of human c-erbB2.
Isotype:	IgG
Format:	This antibody is supplied as tissue culture supernatant containing sodium azide as preservative.
Specificity:	This antibody recognizes a c-erbB-2 protein which is a receptor tyrosine kinase of the c-erbB family. It is closely related in structure to the epidermal growth factor receptor. c-erbB-2 oncoprotein is detectable in a proportion of breast and other adenocarcinomas, as well as transitional cell carcinomas. This antibody crossreacts with Human.
Positive Control:	Breast carcinoma
Cellular Localization:	Cell membrane
Storage:	Store at 2-8°C. Do not use beyond the expiration date stated on the label.

### ASR: Analyte Specific Reagent

DBS will not be held responsible for patent infringement or other violation that may occur with the use of our product

**DBS**

1020 Serpentine Lane, # 114, Pleasanton, CA 94566 Tel: 925 484 3350, Fax: 925 484 3390

Website: [www.dbiosys.com](http://www.dbiosys.com) e-mail: [customersupport@dbiosys.com](mailto:customersupport@dbiosys.com)