

Rabbit Monoclonal Antibody to Human c-erbB-2 Protein

Catalog No.:	RMPD 008
Applications:	Analyte Specific Reagent. Analytical and Performance characteristics are not established. Optimal working conditions must be determined by the end user.
Immunogen:	Recombinant protein encoding extracellular domain of human c-erbB-2.
Clone:	SP3
Isotype:	IgG
Positive Control:	Breast carcinoma
Cellular Localization:	Cell membrane
Storage:	Store at 2-8°C. Do not use beyond the expiration date stated on the label.

ASR: Analyte Specific Reagent

DBS will not be held responsible for patent infringement or other violation that may occur with the use of our product

DBS

1020 Serpentine Lane, # 114, Pleasanton, CA 94566 Tel: 925 484 3350, Fax: 925 484 3390

Website: www.dbiosys.com e-mail: customersupport@dbiosys.com